

HAFILLER®

Beauty Makes You Healthier



Zhejiang Jingjia Medical Technology Co.,Ltd.

Address: 206#, Zhenzhong Road, Xihu District 310030, Hangzhou, China

Tel : 86-571-89987901

Fax : 86-571-88883312

E-mail : biz@hafiller.com

Web: <http://www.hafiller.com>

<http://www.hafiller.com>

SPECIALISTS IN THE DEVELOPMENT, PRODUCTION AND DISTRIBUTION OF Biomedical Materials

HAFILLER®
Beauty Makes You Healthier




ABOUT US JINGJIA MEDICAL

We was founded in 1995, as one of the earliest domestic companies to produce medical sodium hyaluronate gel, is a collection of bio-medical materials research and development, manufacture and sale of state-level high-tech enterprise, medical devices in class III production enterprise.



INTERNATIONAL PRESENCE Distribution networks in



OUR EXCLUSIVE TECHNOLOGIES

WHAT IS HAFILLER

HAFILLER[®] is Medical Cross-linked Sodium Hyaluronate Gel, it is one kind of Soft Tissue Filler and Dermal Filler.

CROSS-LINKED TECHNOLOGY

SAX-HA[™] technology is a new cross linking technology, it can let synthetic hydrogel with hydrophobic on the polymer chains tends to have the strong self- association functionality via the strong hydrophobic association.

This self-association is dynamic.

The injected gel tends to stay and aggregate together to hold the shape and reduce the gel migration in vivo.

Patent No.: ZL 2012 1 0372786.0

PCT No.: PCT/CN2013/080217

ADVANTAGE

No animal origin


No skin test needed

Fully biodegradable and without any rejection


Effect is natural

HAFILLER[®]
Beauty Makes You Healthier



HAFILLER[®] Global Clinical Trial Data Statistics



Effect Maintainace Rate comparison



With a number of clinical trials and doctor's feedback, it showed that Aquaderm cross linked dermal filler have been recognized by global customers. Besides, we will provide continuous clinical safety data support for global patients through a large number of clinical follow-up every year.



Extrusion Force

HAFILLER[®] is extrusion force curve is more smooth than swedish brand, doctors will have a better injection experience with HAFILLER, easier and more smooth.

The frequency distribution of the injection reaction duration during 2 weeks after injected Hafiller

	Total of injection reaction	one day	2-7 days	more than 7 days	more than 14 days
Pain	211 (98.6)	139 (65.0)	72 (33.6)	0 (0.0)	0 (0.0)
Swelling	211 (98.6)	25 (11.7)	186 (86.9)	0 (0.0)	0 (0.0)
Allergy	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)
Congestion	11 (5.1)	6 (2.8)	5 (2.3)	0 (0.0)	0 (0.0)
Induration	102 (47.7)	9 (4.2)	93 (43.5)	0 (0.0)	0 (0.0)

The frequency distribution of the injection reaction duration during 2 weeks after injected other market brand





	Total of injection reaction	one day	2-7 days	more than 7 days	more than 14 days
Congestion	72 (52.2)	7 (5.1)	56 (40.6)	6 (4.4)	3 (2.2)
Rubedo	117 (54.8)	19 (13.8)	68 (49.3)	18 (13.0)	12 (8.7)
Swelling	120 (87.0)	16 (11.6)	84 (60.9)	16 (11.6)	4 (2.9)
Pain	79 (57.2)	29 (21.0)	48 (34.8)	2 (1.4)	0 (0.0)
Tenderness	107 (77.58)	21 (15.2)	78 (56.5)	6 (4.4)	2 (1.4)
Pruritus	42 (30.4)	11 (8.0)	25 (18.1)	6 (4.4)	0 (0.0)
Other	34 (24.6)	7 (5.1)	23 (16.7)	3 (2.2)	1 (0.7)

Enzymolysis in vitro

HAFILLER[®] is enzymolysis in vitro time is less than swedish brand, this date shows that HAFILLER, has longer lasting results.

DERMAL FILLERS



PRODUCT RANGE OVERVIEW

Hafiller Model	Fine Lines	Derm	Derm Deep	Derm Plus	Sub Skin
					
Indication	Fine Lines around eye and mouth. Such as crow's feet and perioral lines.	Moderate wrinkles and skin folds. Such as frown lines. Lip enhancement.	Deep Wrinkles and folds. Such as nasolabial fold. Lip enhancement and face filling.	Breast and buttock augmentation.	Correction of facial contours include cheek, chin and nose.
Injection site					
Particle size	0.10-0.15mm	0.15-0.28mm	0.28-0.50mm	0.50-1.25mm	1.25-2.0mm
Concentration	20mg/ml	20mg/ml	20mg/ml	20mg/ml	20mg/ml
Needle size	30G	30G	27G or 26G	23G or 25G	23G or 25G
Duration	6-9 months	6-12 months	6-12 months	9-18 months	9-18 months

APPLICABLE PARTS



DETAILS
FINE LINES



PRODUCT INDICATIONS

To correct very thin wrinkles, eye around wrinkles and thin wrinkle in facial wrinkles, "crow's feet", wrinkles around the contour of mouth, etc.



Concentration
20mg/ml




Needle
30G



Expected duration
6 to 9 months

1.0ml

2.0ml



PRODUCT INDICATIONS

Correction of medium wrinkles,
skin folds and the red border of the lips.



Concentration
20mg/ml



Needle
30G






Expected duration
6 to 12 months

1.0ml

2.0ml



DETAILS DERM DEEP



PRODUCT INDICATIONS

Medium and deep wrinkles filling,
lips and cheek volume expansion.



Concentration
20mg/ml



Needle
27G or 26G






Expected duration
6 to 12 months

0.5ml

1.0ml

2.0ml

DETAILS DERM PLUS



PRODUCT INDICATIONS

Derm Plus is designed to fill deep wrinkles, recover and form volume of cheeks, cheekbones and chin.



Concentration
20mg/ml



Needle
23G or 25G






Expected duration
9 to 18 months

1.5ml

10ml

20ml

DETAILS SUB SKIN



PRODUCT INDICATIONS

Sub Skin is designed to form volume, correct facial contours, expand volume of cheeks, chin, restore and change of nasal contour, is used to compensate for the loss of volume during skin aging.



Concentration
20mg/ml



Needle
23G or 25G



Expected duration
9 to 18 months

2.0ml

10ml


20ml

BEFORE
AFTER

HAFILLER® FINE LINES: TEMPLE AREA, CROW'S FEET



HAFILLER® FINE LINES: PERIORAL AREA, LABIAL COMMISSURE



HAFILLER® DERM: FORHEAD, GLABELLA



BEFORE
AFTER

HAFILLER® DERM DEEP: CHEEKBONES, NASOLABIAL FOLDS, LIPS



HAFILLER® DERM PIUS: BREAST, BUTTOCK



HAFILLER® SUB SKIN: BREAST, BUTTOCK



DETAILS KNEE INJECTION



PRODUCT INDICATIONS

Application scope: This product assists to prevent or reduce adhesion after abdomen(pelvic) cavity operation.

Storage: Being stored in a place away from light at ambient temperature.



- « Concentration
10mg/ml
- Model
1ml / 1.5ml / 2ml / 2.5ml
3ml / 4ml / 5ml
- © Validity period
2 years

DETAILS OPHTHALMOLOGY



- **Moistening and lubrication**
- **Wound healing**
- **Biocompatibility**

Action:

Maintain normal anterior chamber, protect intraocular tissue, assist pupil dilation, prevent vitreous prolapse, help hemostasis and prevent postoperative adhesion.



« Concentration
14mg / 15mg / 18mg
20mg / 25mg / 30mg

Model
1.5ml

© Validity period
2 years